



# FURIX

FUVCT 900



## Veterinary Fuvision CT Scan



## INTRODUCTION

Vetease 3 has four modes of imaging: CT / DR / Dynamic DR / Dental Panoramic, and matches all the X-ray diagnostic needs of animals.

### 01 600 slices for CT scan

- Super 3D reconstruction

### 02 FOV is 40cm

- Scanning time is 12s

### 02 0.14mm slice thickness

- Big diameter of 800 mm

✓ Unique CT-AIR Technology

✓ Low dosage of radiation, which is just 1/100 of conventional CT scan



✓ 43×43cm large-sized DR FOV



## Clinical Cases



▲ DR



#### Parameters

Function	CT/DR/Dynamic DR/Dental Panoramic/DSA(Upgraded)	Max. CT FOV Width	Φ40cm
Structure	Ring-shaped	Max. CT FOV Length	72cm
Field of DR/Dynamic DR View	43cm × 43cm	Spatial Resolution of CT	1.5 lp/mm
Spatial Resolution of DR/Dynamic DR	25 Strips/vial	Fastest Scanning Mode	12s
Machine Size	2.0m × 1.7m × 1.6m	Gantry Dimension	≥80cm

**Shenzhen Falcon and Furix Lifesciences Co. Ltd**  
**Furix Trade Venture Limited (HK)**

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,  
 Mob:+852 9579 0472 E-mail: sales@myfurix.com

# FURIX