



# FURIX

FDC-201

**Veterinary 3 IN 1 Dry  
Chemistry**



### One device for multiple uses

capable of simultaneously monitoring lipid levels, uric acid, blood glucose, and hemoglobin

### High-definition color screen

4.3inch color screen, clearly display test results

### Accurate results

Stable results and good repeatability

### compact and portable

easy to carry in your pocket , available anytime, anywhere

### Software support

Configure chronic disease management system  
Efficient management of patient groups



multiple uses



Accurate results



color screen



LONG-LASTING POWER

## PRODUCT FEATURE

	Blood Sugar	Uric Acid	Total Cholesterol
measuring range	1.1-33.3 mmol/L	179-1190 $\mu$ mol/L	2.59-12.93 mmol/L
measuring time	5S	20s	3 min
sample type	peripheral whole blood		

	HDL Cholesterol	Triglycerides	Hemoglobin
measuring range	0.39-2.59 mmol/L	0.51-7.34mmol/L	45-256 g
measuring time	3 min	3 min	30s
sample type	Venous or peripheral whole blood		



**Shenzhen Falcon and Furix Lifesciences Co. Ltd**  
**Furix Trade Venture Limited (HK)**

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,  
Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com  
Tel:+852 9579 0472

# FURIX