



FURIX

FD-610C & FD-610C PRO

Furix Dental Unit

High-quality images to assist doctors in accurate diagnosis

APS CMOS sensor enables images with higher signal-to-noise ratio and higher resolution, resulting in higher quality images



3 Way-syringe Handpiece



Low Speed Handpiece



High Speed Handpiece



Ultrasonic Scaler



Suction Handpiece-

Furix Deantal Unit FD-610C

| | |
|-------------------|-------------------|
| Air Flow | 70KL/Min |
| Preddure | 0.8Mpa±0.05mpa |
| Air Tank Capacity | 10L,oil-free,550W |
| noise | 40-55dB |
| N.Weight | 56kgs |
| G.Weigh | 66kgs |

| | |
|--------------|------------------------|
| Product Name | Veterinary Dental Unit |
| Model | FD-610C |
| Voltage | AC220V/50H |
| Current | 3.7 A |
| Power | 550W |
| Speed | 1390/min |



STANDARD CONFIGURATION

Suction HAndpiece * 1 PC
3 Way Syring Handpiece * 1 PC
LED Curing Light * 1 PC
High Speed Handpiece * 1 PC
Low Speed Handpiece * 1 PC
Ultrasonic Scaler 1 PC

Why Choose Our Animal's Dental Care?

Opt for our veterinary dental care featuring the Veterinary Dental Unit FD-610C And FD-610C PRO, renowned for its high-resolution X-ray imaging. This advanced technology ensures accurate diagnosis and effective treatment of your pet's dental conditions, emphasizing their well-being through meticulous care and innovation in veterinary dentistry. Trust us to elevate your pet's dental health with precision and compassion.

Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
Tel:+852 9579 0472 www.myfurix.com

FURIX



Furix Deantal Unit FD-610C Pro

Technical parameters

Processor: i5 three-generation, four-threaded processor
 Hard Drive: 128G high speed solid state hard drive
 Memory: 8G-DDR3
 Screen resolution: 1920 x 1080
 Operating system: Windows 10
 Built-in WIFI, Built-in Bluetooth
 Standard high-definition SUB3.0K oral
 Observation instrument
 Built-in dental observation equipment medical
 Record special software

Display Terminal Board

- ① Power button
- ② USB port
- ③ USB port
- ④ RFIN port
- ⑤ HDMI port
- ⑥ VGA port
- ⑦ USB port
- ⑧ Wireless network port
- ⑨ SPKport
- ⑩ MIC port

Product equipment

| S/N | Name | Qty |
|-----|--------------------------|-----|
| 1 | All-in-one computer | 1 |
| 2 | HD handle (with line) | 1 |
| 3 | DC 12V power adapter | 1 |
| 4 | Power Cord | 1 |
| 5 | HandShank induction Rack | 1 |

| S/N | Name | Qty |
|-----|---------------------------------|-----|
| 6 | Special screw for rack | 8 |
| 7 | Instruction manual | 1 |
| 8 | Simple metal bracket (optional) | 1 |



FD-610C PRO

- 3000k-6000k Linear Adjustable Color Temperature
- 8000-30000lux Adjustable Illuminance, Soft Light Without Hurting Ever
- Infrared Sensor Switch, Power Saving And Convenient
- Arbitrary Angle Adjustment
- Anti-curing Mode, Blue Light Accounts For Orily About 8%

| Features | Specification |
|---------------|---------------------------|
| Input voltage | AC/DC12V/24V |
| Power | 6W(MAX 6W*2) |
| Frequency | 50/60Hz |
| Illumination | 8000-30000Lux(adjustable) |
| Temperature | 3000-6000K(adjustable) |
| CRI | 95 |
| Spot size | 7*140MM (distance 70CM) |

Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
 Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
 Tel:+852 9579 0472 www.myfurix.com

FURIX

Digital oral endoscopic system

Key introduction

Key 1:

- 1 Single - click to photo
long - click to video



- ◎ Patient management: Can help doctors systematize the management of patient files, key words can query patient information.
- ◎ Image capture: The main function is "photo" and "video", the other auxiliary functions are zoom/zoom /actual size/ spread/vertical image/horizontal image/full-screen/video settings, photos and videos taken are automatically saved to a list of images in the lower left corner for comparison and editing.
- ◎ Image contrast: In the image contrast state can also open the camera for real-time shooting, four pictures for comparison, comparison photos can be saved directly to the computer or click on the two-dimensional code for transmission.
- ◎ Image editing: The image can be annotated to facilitate doctors to explain the patient analysis.
- ◎ Report printing: Edit case book, add patient information and real-time examination picture, fill in examination result and doctor's suggestion, print or transfer.
- ◎ System settings: You can choose chinese or english, you can set the name of the hospital and medical records of this report to modify the title of the report icon



Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
Tel:+852 9579 0472 www.myfurix.com

FURIX