




FURIX

Fuvision DRF-32KW

**Fuvision Dynamic Static
DRF Machine.**

Precision Imaging Decoding Every Detail For Life



✓ Patents & Certificates

✓ Superior Performance

✓ Military-Grade Quality

✓ Industry Empowerment



High precision process

Extra thick steel plate
Medical grade spray coating



Sealed design

Dust and hair protection
Industrial-grade PC hardware



Unique Three-door Design

More protective
Easy to open the inside of
the machine for maintenance




Precision, safety and efficiency

Industrial capacitive touch screen
Bedside can be operated by one
person




fVet X-ray Software


Precision, Efficiency, Flexibility




Step1



Step2



Step3



Step4


Shenzhen Falcon and Furix Lifesciences Co. Ltd

Furix Trade Venture Limited (HK)

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,

Mob:+852 9579 0472 E-mail: sales@myfurix.com

FURIX



Specification

X-RAY GENERATOR SPECIFICATIONS

Constant Potential kW ratings	32 kW
kVp Range	40 - 125 kVp
mAs Range	0.1mAs~125mAs
mA Range	Photography mode: 10~400 mA, Continuous fluoroscopy mode 0.5~10mA Pulse fluoroscopy mode: 0.5~100mA
Exposure Time Range	0.001-20 s

Tube

X-Ray tube	Rotating Anode / (0.6 / 1.2 mm)
Heat capacity:	300KHU
Anode input power:	Small focus 20 to 22 KW/Large focus 40~42 KW

Rack

Table - top type	Floating - 4 way table-top
Table top lenght	Standard 140 cm, Optional: 120cm
Maximum load patient weight	200 Kg
Brakes	Electromagnetic

Flat panel detector

Scintillator fluorescent material	CsI
Pixel spacing	139µm
Maximum load patient weight	200 Kg
Protection class	IP65
Pixel matrix	3072 × 3072 pixels
A/D Conversion	16 bits

X-rays that can capture animal dynamics

Iteration from pictures to videos

Observe the animal's condition more clearly in its dynamic state and make a more accurate diagnosis

Voice-controlled exposure

Capture the images you need efficiently and accurately

One-key Power On/Off

Fast power on/off
Streamline startup operationsto avoid increased errors

Shenzhen Falcon and Furix Lifesciences Co. Ltd
Furix Trade Venture Limited (HK)

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob:+852 9579 0472 E-mail: sales@myfurix.com

FURIX