



FURIX

FOC-5L

Furix Oxygen Concentrator

DAILY OXYGEN PROTECTION AND HEALTH CARE



FOC-5L



Pet Healthcare Oxygen Concentrator

01

High flow



1-5L/min adjustable flow rate, guarding your pet's oxygen health.

02

Oxygen purity



90%~96% high concentration oxygen, the oxygen outlet is cleaner and more stable.

03

Multi-function



Oxygen generation and atomization, stable and continuous oxygen supply.

04

Low noise



Multiple noise reduction without disturbance, breathing oxygen at ease is super comfortable.

05

Intelligence



One-button oxygen production is easy to operate, easy to use voice broadcast.

Product Specification

Product application

For pets health care
Contraindications: Oxygen poisoning, oxygen allergy patients are prohibited.

Principles of oxygen production

This machine takes ~ 220V power supply as power source, air as raw material, adopts high quality molecular sieve, and produces high purity oxygen by pressure swing adsorption separation method (PSA method) at room temperature.

Security system

- Current overload or connection loosening, shutdown;
- High temperature of compressor, shutdown;
- Low oxygen concentration alarm;

Normal operating environment

Ambient temperature range: 10°C ~ 40°C;
Relative humidity range: ≤ 80%;
Atmospheric pressure range: 860 hPa ~ 1060 hPa.
Note: When the storage environment is lower than 5°C, the equipment shall be placed at normal working temperature for more than four hours before use.

Storage and transportation environment

Ambient temperature range: -40°C ~ +70°C;
Relative humidity range: 10% ~ 100%, without condensation
Atmospheric pressure range: 860 hPa ~ 1060 hPa.
Note: The oxygenator shall be stored in a well-ventilated room where there is no strong sunlight, no corrosive gas, and it shall be transported to avoid sharp shock and inverted recumbency.

FEATURES

- The oxygenator is composed of host, flow indicator and humidification cup. (Configure the atomizing device with medical device registration certificate.)
- Fully plastic enclosure, safe and reliable.
- It has real voice function and more intuitive operation.
- Color screen display, clearer.
- Remote control function: The machine can be operated remotely and is easy to use.
- With accumulation function: display the total working time through the display screen.
- It has the function of timing shutdown, which is easy to use.
- There is power failure alarm function.
- It has the function of releasing oxygen-rich ions.
- The compressor is equipped with a thermal protector to better ensure the safety of the compressor and complete machine.
- With flow regulation function:
Maximum Recommend Flow Rate: 3L/min 5L/min
- It has atomization function.
(Applicable to models with atomization function)
- It has oxygen concentration monitoring alarm function.

Dimension

Overall dimension: 36×30.5×58.5 (cm)

Weight

Net weight: 17kg

Main technical indicators

Oxygen concentration: 93% ± 3%, > 5L: 30% - 90%

Package including

- Oxygen concentrator: 1 set
- Fuse: 1 pc
- User manual: 1 pc
- Oxygen inhalation tube: 1 pc
- Oxygen mask: 1 pc
- Remote controller: 1 pc
- Silicone tube: 1 pc
- Power cord: 1 pc



Imported molecular sieve gives you a full sense of security with natural enjoyment!

Furix

— Furix Oxygen Concentrator FOC-5L



Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,
Mob: +852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com
Tel: +852 9579 0472 www.myfurix.com

FURIX