



# FURIX

FOC series

Oxygen Generator System



# Oxygen Generator System

Furix innovative oxygen generating system is designed for hospitals and clinics with 5-100 beds. The system includes PSA oxygen generator, buffer tank, booster pump, high pressure oxygen compressor, etc. It represents a higher medical standard and true sense of innovation.

## 01 Features

- It can be installed anywhere without site requirements.
- Low noise, it can be used in ICU and operation rooms.
- It offers an economic, reliable and safe alternative to traditional oxygen system.
- The design is made for 24/7 operation.
- It is equipped with automatic start&stop function, enabling the system to start and stop automatically.

## 01 Application



- For hospital with 5-100 beds.
- For anaesthesia machine and ventilator in operating room and ICU.
- For high pressure cylinder filling plant.



## Specification



Model	FOC-15	FOC-20	FOC-30	FOC-40	FOC-60
FLOW RATE	1-15L	1-20L	1-30L	1-40L	1-60L
SOUND LEVEL	43DB	48DB	58DB	58DB	60DB
POWER CONSUMPTION	600W	1200W	1650W	2200W	3300W
NET WEIGHT	90 KGS	95 KGS	140 KGS	160 KGS	245 KGS
DIMENSION	530*440*780MM	530*440*1380MM	760*540*1380MM	760*540*1380MM	760*540*1980MM
OXYGEN CONCENTRATION	93% ( ±3% )				
OUTLET PRESSURE	6-10 Psi/ 20 Psi/ 60 Psi				
ALARMS	Power failure alarm;high/low pressure alarm;low oxygen concentration alarm				



**Shenzhen Falcon and Furix Lifesciences Co. Ltd**  
**Furix Trade Venture Limited (HK)**

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,  
Mob:+852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com  
Tel:+852 9579 0472

# FURIX