



# FURIX

FN1912

## Veterinary Circulation Water Warming System

## Touch Operation Convenient and Simple



One-button switch Off/On




Sleeping Mode



Mites Killing



Water shortage reminder



Constant temperature setting, warm water circulation. Keep moist and comfortable without drying.

### Veterinary Water Warming System

Rated Voltage	AC110-120/AC220-240V
Rated Frequency	50/60Hz
Rated Power	280W
Temperature Range	Normal: 10 - 65°F/ Mite killing: 140°F (60°C)
Host Size	19*19*17cm
Mattress Pad Size	Mattress pad size (100*55cm)
Water Tank Capacity	800mL
Host db	Around 19db



### Standard Package

- ① 1 Main Host
- ① 1 Queen Size Anti-scratch Mattress pad size (55\*55cm)
- ① 1 King Size Anti-scratch Mattress Pad (55\*100cm)



Hydrating  
Nourishing



Heating  
Insulation



Mites  
Killing



Precise  
Temperature  
Control



Rapid  
Heating



Protection  
Safety

## Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG,  
Mob: +852 9579 0472 E-mail: vet.furix@gmail.com, sales@myfurix.com  
Tel: +852 9579 0472 www.myfurix.com

# FURIX