




FDX-1  
Pet Dental Portable  
X-ray Machine

FURIX


Product specification

Tube voltage	60KV	Battery	2200mAh
Tube current	2.0mA	Environmental temperature	5°C-40°C
Exposure time	0.1-1.2s	Relative humidity	<75%RH
Frequency	30KHz	Net weight	1.9KG
Rated power	120VA	Gross weight	4.2KG
Distance from focus to skin	118mm	Dimensions (mm)	168*158*130
Tube focus point	0.8mm*0. 8mm	Package size (mm)	365*330*275


Product Features




Small size, strong penetration



Easy to operate with just one click



Child lock function safety protection




Built in lithium battery

Product Features

- 01 0.4mm focusing technology provides clear and high-quality images to meet various diagnostic needs
- 02 The entire machine adopts advanced technology of imported high-frequency tubes to ensure the stability and safety of ray output
- 03 High magnification power battery, constant charge resistance, fully charged, sustainable for 200 shots
- 03 Pet specific dental profession, making you more relaxed and professional

Product Highlights

- SLR design with simple filming
- Short exposure time
- Take and take photos anytime, easy to operate
- Strong ray penetration
- New pet for filming: exclusive pet parameter design
- One click settings for filming



Furix Trade Venture Limited

Address: 6/F MEGA CUBE, WANG KWONG ROAD, KOWLOON, HONG KONG  
E-mail:vet.furix@gmail.com, sales@myfurix.com Web:www.myfurix.com  
Tel:+852 9579 0472 Mob:+852 9579 0472

FURIX